


# Salzburger Landesklinikum - St. Johannis Spital

- Standard-Monovetten für Zentrallabor -

## Blutentnahmesystem - S-Monovette® für Zentrallabor


SAP-Nr.: 70023287



Li-Hep-Gel

Plasma-Gel  
4,7 ml


SAP-Nr.: 70009395



Serum

Serum  
5,5 ml


SAP-Nr.: 70009379



Coagulation

Citrat 1:10  
3,0 ml


SAP-Nr.: 70009374



EDTA K

EDTA K  
2,7 ml


SAP-Nr.: 70009383



Glucose

Fluorid  
2,7 ml


SAP-Nr.: 70009391



BSG / ESR

BSG  
2,0 ml


SAP-Nr.: 70009398



Citrate/Buffer 9 NC/PFA

PFA  
3,8 ml


SAP-Nr.: 70009399



Blood gas

Blutgas  
2,0 ml

Vom Labor zu beziehen!



ThromboExact / 2.7 ml  
Pseudothrombocytopenie

ThromboExact  
2,7 ml

Blutentnahme-Röhrchen von anderen Firmen hier nicht dargestellt